ML-0051

PHANTOM LED LIGHT SOURCE

Operator Manual

TeDan Surgical Innovations
12320 Cardinal Meadow Drive, Suite 150
Sugar Land, TX 77478
Phone: 713-726-0886
Fax: 713-726-0846

Sunoptic Technologies®
6018 Bowdendale Avenue
Jacksonville, FL32216 USA

AJW Technology Consulting GmbH
Breite Straße 3
40213 Düsseldorf, Germany

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Table of Contents

Symbology	3
Introduction	4
Intended Use	4
Contraindications	4
Patient Target Group and Intended User(s)	4
Adverse Events	5
Clinical Benefits	5
Adverse Effects	5
Prior to Use	5
Expected Life of Device	5
Storage Conditions	5
Function of Design	6
Warnings / Cautions	6, 7
Warranty	7
Specifications	8
Certifications	9
Device Overview	10, 11, 12, 13, 14
Device Set-Up and Operation	15
Cleaning and Disinfection	16
Maintenance, Service/Repair	16
Troubleshooting	17
End of Product Life	17

1. Symbology

Symbol	Description	Symbol	Description
 or  YYYY-MM-DD	Manufacture Manufacture w/Date of Manufacture NOTE: Date may also be placed on right side of factory symbol.		Date of Manufacture (YYYY-MM-DD)
	“Authorized Representative”, European Community		“Authorized Representative”, Switzerland
	Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner		European CE Mark
	Instructions for Use		Medical Device
	Supplied Non-Sterile		Unique Device Identifier
	Note		Caution: Hot Surface
	Caution		Caution: High Voltage
	Warning		Keep Dry
	Do Not Use if Damaged		Special Electrical / Electronic Disposal required.
	Storage / Shipping Temperature		Storage / Shipping Humidity
	Barometric Pressure		Product Safety Tested CSA Canada, US
	Alternating Current		Equipotentiality
	Type BF		Protective Earth (Ground)
	Power “ON” Symbol		Power “OFF” Symbol

2. Introduction

Congratulations on the purchase of your new PHANTOM ML-0051 LED Light Source!

This user-friendly LED light source is a high efficiency illumination system utilizing state-of-the-art superior LED illumination technology. It offers a variety of features such as:

- 6500 K daylight brightness for good color definition
- Quiet operation
- Long life, 50,000 hrs.
- Compact and light weight
- Turret which adapts to various types of light guides
- Pulse-width Modulation (PWM) Electrical Dimming

In short, you have chosen the best and we would like to make sure you receive the optimal results with your new LED Light Source by using it correctly.

This Operator Manual will help you to install the device and optimally integrate it with other components of your system. It will also instruct you how to operate the LED Light Source and how to keep it clean. It will give you maintenance and service guidelines as well as recommendations for best performance results.

	<p>This light source is designed to work with TeDan Surgical Innovations the proprietary TSI lights cables light cable only to be used during spine, hip, neuro, cardiac, vascular and thoracic surgeries.</p> <p>Contact TeDan Surgical Innovations for additional information.</p>
---	--

	<p>No special intervention is necessary in the event of device failure. As such, this device has no essential performance, as defined by IEC 60601-1.</p> <p>Should this device fail to operate, a suitable backup should be available for any procedure where it may be used.</p>
--	--

3. Intended Use

The LED light source is intended to provide illumination to a surgical site during minimally invasive procedures including, but not limited to, arthroscopy (orthopedic surgery), laparoscopy (general and gynecological surgery), and endoscopy (general, gastroenterological, and ENT surgery).

Light is typically transmitted from the illuminator through a fiber-optic cable and typically to a scope or other light guide.

	<p>This LED light source is intended to be used in a controlled surgical setting environment by qualified medical personnel.</p> <p>The LED light source is provided in a non-sterile state and is not intended to be sterilized. Refer to Cleaning Instructions provided in this manual.</p>
---	---

4. Contraindications

This device is not intended to be used for monitoring, diagnostic, or other life support functions. It is not intended to compensate for injury, handicap, replacement or modification of anatomy, or control of conception.

	WARNING	<p>This device should <i>NEVER</i> be used in ocular surgery or in a surgical procedure requiring direct illumination of the eye.</p>
---	----------------	--

5. Patient Target Group and Intended User(s)

This device is not intended for a particular target group or groups, as the device is intended to be used in a surgical setting for illumination of a surgical site under the supervision and direction of a physician only.

6. Adverse Events / Serious Incident – Medical Device Reporting

Per the FDA, an Adverse Event or Serious Incident is defined as:

Caused or contributed means that a death or serious injury was or may have been attributed to a medical device, or that a medical device was or may have been a factor in a death or serious injury, including events occurring as a result of:

- Failure
- Malfunction
- Improper or inadequate design,
- Manufacture
- Labeling, or
- User error.

	In the event of an Adverse Event or Serious Incident, the end-user should contact the distributor immediately.
--	--

7. Clinical Benefit - Definitions

Regulation (EU) 2017/745

Article 2(53) MDR defines clinical benefit as the positive impact of a device on the health of an individual, expressed in the terms of a meaningful, measurable, patient-relevant clinical outcome(s), including outcome(s) related to diagnosis, or a positive impact on patient management or public health.

FDA

Clinical benefit is a favorable effect on a meaningful aspect of how a patient feels (e.g., symptom relief), functions (e.g., improved mobility) or survives as a result of treatment. • Clinical benefit may be measured as an improvement or delay in the progression of a disease or condition (as manifested by how a patient feels/functions).

	The manufacture of this device makes no claims as to any Clinical Benefit to be gained by the use of this product.
--	--

8. Adverse Effects

The manufacture is not aware of any Adverse Effects experienced from the use of this device when properly utilized per the devices intended use. Refer to Section 6 of this document for additional details.

9. Prior to Use

It is recommended that the end-user perform a functional check of the light source to ensure proper operation prior to use of a surgical procedure.

10. Expected Life of Device

This device has an expected service life of 3-years.

11. Storage Conditions

- Temperature Range: -4 to 140°F (-20 to 60°C)
- Humidity Range: 30 to 95% RH non-condensing
- Barometric Pressure: 700 to 1060 hPa

12. Function(s) of Design

The LED Light Source is a standalone illuminator that produces 1000lm (typical) in the visible spectrum range.

The LED light source is typically used with a fiberoptic light guide that connects to the light port located on the front panel.

Light intensity output is controlled via rotating dimming knob located on the front panel of the console.

13. Warnings / Cautions

Use of this equipment may present hazards to the user and/or patient. Before operating this device, please read this operating manual thoroughly and follow all warnings, cautions, and instructions for use.

	WARNINGS	WARNING symbol indicates conditions that may cause risk to the safety of the patient or user. Failure to adhere to the follow warnings may result in injury to the patient or end user.
---	-----------------	---

- The illuminator produces highly concentrated light. Avoid shining light beam into eyes or looking directly into the light beams at the ends of connected instruments and/or light guides. To reduce risk, always install a fiberoptic cable when the illuminator is OFF or when the dimming knob is turned to the fully dim position. When not using the device, it is advised to fully dim the illuminator.
- Qualified personnel must determine a safe working distance between the ends of connected instruments and/or light guides and the patient for each application. There is a risk of patient injury if a light guide or instrument connected to the light source comes to close to the patient.
- User is responsible for determining if interruption of light output will create an unacceptable risk. Having a backup illuminator is advised.
- User is responsible for providing backup lighting systems for your application when using this device.
- Not suitable for use in presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- For endoscopic procedures: the illuminator should only be used with type BF endoscopic instruments which have been certified to IEC 60601-1 and IEC 60101-2-18.

This symbol indicates type BF equipment.

- All devices and/or instruments connecting to the illuminator must be classified as medical equipment. It is the user's responsibility to ensure that all equipment used with this device meets all applicable standards such as IEC 60601-1.
- To prevent fire and/or electric shock, do not open or expose the illuminator to liquids.
- The fiber optic light guide connected to the LED light source may become hot during use. Allow adequate time to for end tip to cool before removal from illuminator.
- Instruments and/or light guides connected to the illuminator must be **Non-Conductive**. There should be no conductive shielding or any conductive connection between the illuminator and the patient. Such connections present a risk to patient safety.
- Instruments and/or light guides should be clean and dry before being connected to the illuminator.
- **DO NOT** modify the equipment without authorization from the manufacturer.
- The illuminator is provided non-sterile and is not intended to be sterilized.
- Use only the power cord supplied with illuminator or medically approved power cords with less than 200mΩ of ground impedance and less than 16ft (<5m) of length. If unauthorized cables are used, the device may have increased electromagnetic emissions and/or decreased electromagnetic immunity which may result in improper operation.

- This device meets CISPR 11 Class A limits and is suitable for use in a hospital and industrial environments. If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
- Performance of this device may be affected in proximity of another device and/or equipment capable of producing high levels of RF emissions. The device should be used no closer than 12 inches (30 cm) to any part of RF equipment including cables. In the event performance of this device is affected due to high levels of RF emissions, relocation of the suspected device and/or equipment producing high levels of RF emissions, or the headlight system may reduce or eliminate the problem.

	CAUTIONS	Indicates risks of improper use and/or damage to the equipment. Failure to adhere to the following cautions may result in loss of function or product damage.
---	-----------------	---

- Before each procedure, carefully check the illuminator for damage. ***DO NOT*** use if damaged.
- User should verify the light guide end tip and the active illumination port are of the same type before insertion. ***DO NOT*** attempt to force an end tip into an incorrect port.
- All service and repair must be performed by the manufacturer or qualified service technicians.
- Ensure that the air vents located on the illuminator are not obstructed to allow the device to receive the necessary cooling to prevent an overheating.

	Any serious incident occurring in connection with the device must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or the patient is established.
---	---

14. Warranty

	Upon receipt of the product, the device should be carefully inspected for visible damage. If any defect is discovered, notification must be given immediately to the distributor.
---	---

This LED light source carries a 3-year warranty from the date of shipment on workmanship and all defects of material. Should your product prove to have such defects within three year of the shipment, **TeDan Surgical** will repair or replace the product or component part without charge. Should your LED Light Source product(s) need servicing under this warranty, please contact **TeDan Surgical** for return authorization documentation. You should carefully pack unit in a sturdy carton and ship it to the factory. Please include a note describing the defects, your name, telephone number and a return address. Warranty does not cover equipment subject to misuse, accidental damage, normal wear and tear or if transferred to a new owner without authorization from **TeDan Surgical**. This warranty gives you specific legal rights and you may also have other rights that vary from state to state.

15. Specifications

PARAMETER	VALUE
Light Source Type	LED (Light Emitting Diode)
Color Temperature	6500 K (typ.)
LED life	50,000 hours (typ.)
CRI	70 (typ.)
Brightness Control	PWM (Pulse-width Modulation) 0-100% dimming
Light guide adaptor	Rotating turret with STORZ, ACMI, WOLF, and OLYMPUS (if equipped)
Input Power	100-240 VAC, 50/60 Hz 0.5 – 1A (max.)
Fuses	¼" x 1-1/4" (6.4 x 31.7mm), 250V, 2A, Time Delay
Operating Conditions	32 to 104°F (0 to 40°C), 30 to 85% RH non-condensing, 700 to 1060 hPa
Storage Conditions	-4 to 140°F (-20 to 60°C), 30 to 95% RH non-condensing, 700 to 1060 hPa
Dimensions	4.65" x 3.32" x 7.25" (W x H x D) 11.8 x 8.4 x 18.4 cm (W x H x D)
Weight	2.5 lbs / 1.14 kg

16. Certifications

PARAMETER	VALUE
System Classification	FDA Class I, Device Listing D095692, 510(k) exempt EU Class I, Active device per Annex IX, rule 1
Isolation	Type BF
EMC Certifications	CISPR 11 Class A, IEC 60601-1-2 4 th Edition Electrostatic discharge: ± 8 kV contact, ± 15 kV air Radiated RF EM Fields: 3 V/m, 80 – 2700 MHz EFT / Burst: ± 2 kV, ± 1 kV signal lines, 100 kHz SURGE: ± 0.5 , ± 1 , ± 2 kV unsym Conducted disturbance: 3 V 150 kHz – 80 MHz and 6 V in ISM bands Power frequency magnetic fields: 30 A/m Voltage dips: 0% U_n / 0.5 cycles at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°, 40 % U_n / 5 cycles, 70 % U_n / 25 cycles Voltage interruptions: 0% for 5000 ms Proximity fields: in accordance with EN 60601-1-2: 2015 table 9
CE Marking	Regulation (EU) 2017/745
Other Approvals	CAN/CSA-C22.2 No. 60601-1-08, TC 2:2011 (Corrigendum 2) CAN/CSA-C22.2 No. 60601-2-18:11
Degree of protection against harmful ingress of water	IPX-0; no protection.
Degree of safety in the presence of Flammable Anesthetics	Equipment is NOT suitable for use in the presence of flammable anesthetics.
Mode of operation	Continuous


17. Device Overview

FRONT PANEL

No.	Description	Function
1.	Power switch	Used to turn light source "ON" or "OFF"
2.	Light guide turret	Connection for TSI fiber optic cable
3.	Intensity Control knob	Controls the light output intensity

REAR PANEL

No.	Name	Function
1.	AC Mains Inlet	Accepts AC Power cord
2.	Fuse	2A, 250V (1/4" x 1 1/4") time delay
3.	Grounding Connector	For potential equalization
4.	Venting grids	Allows adequate ventilation for cooling function of the light source


TOP PANEL

Front


No.	Name	Function
1.	Label	TeDan Surgical Innovations Label

LEFT SIDE CHASSIS

No.	Name	Function
1.	Label	Distributor Label Information

BOTTOM PANEL

No.	Name	Function
1.	Product Label	Product Label containing: Model Number, Serial Number, Electrical Ratings, Manufacturer Name, Regulatory Marks and FDA "Rx Only" Symbol, UDI Information and Manufacture Date.
2.	Manufacture Label	Legal Manufacture Information

18. Device Setup and Operation

	WARNINGS	Indicates risks to the safety of the patient or user. Failure to follow warnings may result in injury to the patient or user.
<ul style="list-style-type: none"> ▪ DO NOT use the device in any environment with explosive or flammable gases. ▪ Avoid placing the device in an area where the illuminator may be splashed with liquids. ▪ DO NOT obstruct the exhaust or cooling vents of the illuminator. User must ensure that environment air temperatures surrounding the unit are within the allowable limits. ▪ Ensure that any cables and/or instruments connected to the device are not obtrusive to the procedure. ▪ To prevent electric shock, connect power cords of peripheral equipment through medical isolation transformers. ▪ To avoid adverse electromagnetic effects, DO NOT operate this equipment near high RF energy equipment. 		
	CAUTIONS	Indicates risks of improper use and/or damage to the equipment. Failure to follow cautions may result in loss of function or product damage.
<ul style="list-style-type: none"> ▪ Use only power cords provided with the unit or power cords approved for medical use. 		
	When using a medical isolation transformer, ensure that the transformer has sufficient power ratings. Ensure that the power cord is connected to mains power with a three-prong plug.	

18.1 Device Setup

- Place the Illuminator on a stable surface such as a cart, counter, stand, etc.
- Make sure the power switch is in the “OFF” position.
- Connect AC power cord to the appliance inlet located on the rear panel of the light source.
- Connect a fiberoptic light guide to the turret, ensuring that the fiberoptic end tip type matches the turret port.

18.2 Operation

- After the power cable and light guide are properly connected, rotate the dimming knob to the lowest output position.
- Turn the illuminator “ON” using the power switch located on the front panel. The power switch indicator light, located within the power switch, should illuminate.
- Ensure the connected light guide is safely positioned and adjust the brightness to the desired level by rotating the intensity control knob.

19. Cleaning and Disinfection

	WARNINGS	Indicates risks to the safety of the patient or user. Failure to follow warnings may result in injury to the patient or user.
<ul style="list-style-type: none"> ▪ Ensure that the illuminator is de-energized and disconnected from mains power before attempting to clean and disinfect. ▪ DO NOT use strongly caustic or acidic cleansers such as “Clorox” hypochlorite bleach, ammonia, muriatic acid, or similar products. ▪ DO NOT use acetone, methyl ethyl ketone, or halogenated / chlorinated hydrocarbon solvents or cleansers containing any of these restricted compounds. ▪ The illuminator is not sterilizable. DO NOT attempt to autoclave the device. 		

19.1 Recommended Cleaning Process:

- Use a fresh CaviWipe to clean surfaces for a minimum of 90-seconds.
 - Repeat step above as needed until all device surfaces are visibly clean.
- Wet a low-lint cloth with critical water (DI or RO water) and rinse all surfaces for a minimum of 30-seconds.

	Follow all applicable bloodborne pathogen procedures as required by OSHA and/or your hospital when cleaning and disinfecting the product.
---	---

20. Maintenance & Service/Repair

	WARNINGS	Indicates risks to the safety of the patient or user. Failure to follow warnings may result in injury to the patient or user.
<ul style="list-style-type: none"> ▪ Ensure the light source is turned “OFF” and the device is unplugged from the Mains Power Source (plug removed for power outlet) prior to replace the fuse. 		

20.1 Fuse Replacement

- Turn the light source “**OFF**” and unplug power cord.
- Remove the defective fuse by rotating the fuse holder cover and extracting the fuse from the holder.
- Replace blown fuse with ¼” x 1-1/4” (6.4 x 31.7mm), 250V, 2A, Time Delay fuse.
- Reinstall fuse holder cover.

20.2 Service/Repair

- It is recommended that this device be returned to the distributor for repair purposes. The device should be inspected and evaluated by a qualified technician.
- If the device is no longer in Warranty, an estimate of repair charges will be submitted for approval. Payment must be received before any repair(s) are completed.

	For additional Information or for Repair/Service, contact the distributor per the following information: In the US: 713-726-0886 FAX number: 713-726-0846 e-mail: Sales@tedansurgical.com
---	---

	Any serious incident that has occurred to the end-user or patient, in relation to the device, should be reported to the distributor and/or the competent authority of the member state in which the end-user and/or patient is located.
---	---

21. Troubleshooting

Problem	Solution(s)
The device power switch is not illuminated when in "ON" position.	<ul style="list-style-type: none"> • Check that the AC power cord is properly connected. • Check the unit fuses and replace if necessary. • Contact Distributor.
The device power switch is illuminated but there is no light output from device.	<ul style="list-style-type: none"> • Turn intensity control knob clockwise to increase light output intensity. • Check to ensure fiber optic cable is inserted properly. • Contact Distributor.

22. End of Product Life

For disposal, observe all relevant regulations and/or laws applicable to your location.

- In accordance with the European Waste from Electrical and Electronic Equipment (WEEE) directive, end-users are encouraged to recycle this product whenever possible.
- In the US, perform a Web search for a list of recyclers in your area.

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